Pulsed Azidohomoalanine Labeling in Mammals (PALM) detects changes in Liver-specific LKB1 knockout mice

Daniel B. McClatchy¹, Yuanhui Ma¹, Chao Liu², Benjamin D. Stein¹, Salvador Martínez-Bartolomé¹, Debbie Vasquez³, Kristina Hellberg³, Reuben J. Shaw^{3,4}, and John R. Yates III^{1*}

¹ Department of Chemical Physiology and Molecular and Cellular Molecular Biology, The Scripps Research Institute, La Jolla, CA 92037

²Key Lab of Intelligent Information Processing, Institute of Computing Technology, Chinese Academy of Sciences, Beijing 100190, China

³Molecular and Cell Biology Laboratory, ⁴Howard Hughes Medical Institute, The Salk Institute for Biological Studies, La Jolla, CA 92037

Table of Contents:

Table S-1: List of the amino acids (g/kg) in the standard, PALM, and control diets

Figure S-1. Verification of removal of LKB1 from mouse livers.


Figure S-2. Correlation of other technical replicates of quantitated peptides from NSP not shown in Figure 4B.

All proteins identified and quantified are deposited at http://sealion.scripps.edu/pint/?project=f302f024b4d5aa1f, where the data can be queried and is linked to additional protein information. The data is organized by experiments and figures outline in the main text.

^{*} Corresponding author: jyates@scripps.edu


	DIET			
Amino Acids (g/kg)	Standard		PALM	Control
AHA		0	2	0
Alanine		3.5	9.23	9.23
Arginine	1	2.1	3.7	3.7
Asparagine		6	10.25	10.25
Aspartic Acid		3.5	12.06	12.06
Cystine		3.5	3	3
Glutamic Acid		40	49.45	49.45
Glycine	2	3.3	28.12	28.12
Histidine		4.5	3.5	3.5
Isoleucine		8.2	4.25	4.25
Leucine	1	1.1	7.5	7.5
Lysine		18	5.5	5.5
Methionine		8.2	0	2
Phenylalaine		7.5	7.5	7.5
Proline		3.5	3.5	3.5
Serine		3.5	12.32	12.32
Threonine		8.2	4.25	4.25
Tryptophan		1.8	1.25	1.25
Tyrosine		5	5	5
Valine		8.2	5.5	5.5

Table S-1: List of the amino acids (g/kg) in the standard, PALM, and control diets



Supplementary Figure 1

Figure S-1. Verification of removal of LKB1 from mouse livers. The same mouse livers that were analyzed by PALM analysis were prepared for immunoblot analysis. LKB1 was detected in the WT livers but not the KO livers. Upon deletion of LKB1, loss of AMPK(P-AMPK) phosphorylation on Thr172 was observed as previously reported¹². Total AMPK(T-AMPK) and beta-actin were employed as loading controls. The asterisk denotes the full-length protein.



Supplementary Figure 2

Figure S-2. Correlation of other technical replicates of quantitated peptides from NSP not shown in Figure 4B.