

Parallel accumulation – serial fragmentation (PASEF): Multiplying sequencing speed and sensitivity by synchronized scans in a trapped ion mobility device

Florian Meier^{1‡}, Scarlet Beck^{1‡}, Niklas Grassl¹, Markus Lubeck², Melvin A. Park³, Oliver Raether², Matthias Mann^{1}*

¹ Proteomics and Signal Transduction, Max-Planck-Institute of Biochemistry, Am Klopferspitz 18, 82152 Martinsried, Germany

² Bruker Daltonik GmbH, Fahrenheitstr. 4, 28359 Bremen, Germany

³ Bruker Daltonics Inc., Massachusetts 01821, Billerica, United States

SUPPORTING INFORMATION

- 1) **Supplementary Figure S-1: Annotated PASEF-MS/MS spectra of the peptides shown in Figure 4.**




Figure S-1. Annotated MS/MS spectra of four peptides measured with the PASEF method as shown in Figure 4.