

Capture and printing of fixed stromal cell membranes for bioactive display on PDMS surfaces

Jungwoo Lee¹, Jennifer B. Wang¹, Francesca Bersani², Biju Parekkadan^{1,3*}

Supporting Information



Supporting Figure 1. Morphological comparison fixed-PDMS transferred human dermal fibroblasts, MDA-MB-231 breast cancer cells and DU-145 prostate cancer cells under optical microscope.




Supporting Figure 2. (A) Autofluorescent confocal images of remnant MDA-MB-231 cellular bodies on a tissue culture substrate after cell-PVA film stripping process. (Ex: 458nm / Em: 480-520nm(Green) / Emission: 580-620nm(Red)) **(B)** 3D reconstructed confocal images.


¹ Department of Surgery, Center for Engineering in Medicine and Surgical Services, Massachusetts General Hospital, Harvard Medical School and Shriners Hospitals for Children in Boston, MA, 02114, USA

² Cancer Center, Massachusetts General Hospital, Harvard Medical School, Charlestown MA, 02129, USA

³ Harvard Stem Cell Institute, Boston, MA, 02155, USA



Supporting Figure 3. Nucleus and plasma membrane staining of fixed-PDMS transferred MDA-MB-231 breast cancer and PC-3 prostate cancer cells.



Supporting Figure 4. Fixed-PDMS transferred micro-patterned MDA-MB-231 breast cancer cells.