

Supporting Information

Solution pH that minimizes self-association of three monoclonal antibodies is strongly dependent on ionic strength

Shantanu V. Sule,[†] Jason K. Cheung,[‡] Valentyn Antochshuk,[‡] Amardeep S. Bhalla,[‡] Chakravarthy Narasimhan,[‡] Steven Blaisdell,[±] Mohammed Shameem[‡] and Peter M. Tessier^{‡,}*

[†]Center for Biotechnology & Interdisciplinary Studies, Department of Chemical & Biological Engineering, Rensselaer Polytechnic Institute, Troy, NY 12180.

[‡]Sterile Product Development, Merck Research Laboratories, Summit NJ 07901

[±]Analytical Characterization, Merck Research Laboratories, Union NJ 07083




FIGURE S1. Effect of injection concentration on antibody retention. mAb1 and mAb2 were buffer-exchanged to pH 4.4 and 50 mM NaCl, and injected at concentrations of 3, 5 and 7 mg/mL into columns with immobilized antibody.